

Quick Test™ - Corrosion

Quick Test™ - Corrosion Strips

Fast and accurate testing of presence of Corrosion in hydronic heating systems.

- Easy to use
- Fast and accurate testing
- Onsite testing with no initial set-up requirements
- Results in a matter of seconds
- Cost effective alternative to other solutions on the market
- Available on airtight waterproof containers
- Master packs of 18

Physical and chemical properties

Appearance:	Plastic Strip
Colour:	White
Odour:	N/A

Product uses

Fast and accurate testing of corrosion levels inside hydronic heating systems.

Dissolved iron present in the hydronic heating system's water can indicate a corrosion issue, this can lead to fouling of components and reduction in efficiency.

How to use

1. Drain small sample of system water from the system
2. Dip the Strip in water for 2 -3 seconds
3. Shake off excess liquid
4. Leave for 30 seconds as the water soaks into the pad
5. After 30 seconds compare result with color scale on the side of the pot to indicate levels of corrosion inside the system.

Environmental

ADEY is a carbon neutral company and cares about the environment. The company's award-winning product innovations contribute significantly to reduced carbon emissions, energy and water usage.

Quick Tests are classified as non-hazardous and non-toxic.

Quality

ISO 9001

Packaging, handling and storage

Available in air-tight plastic container.


Keep out of reach of children. Avoid contact with eyes and skin. In case of contact with eyes, bathe the eye with running water for 15 minutes. In case of contact with skin, wash with plenty of soap and water. Store in a cool, well-ventilated area. Dispose of in accordance with local regulations.

DIRECTIONS

Drain small sample of system water and dip a test strip into solution for 2-3 seconds.
Shake off excess liquid.
After 30 seconds, compare the result with the color scale below.



Corrosion Quick Test™



Iron 0-100ppm

PPM (mg/L)

Without MCI®

With MCI®

Fast and accurate testing of corrosion levels

Results within 30 seconds

Contains 25 test strips






Test

adey.com/us

Product ref: CH2-03-06193
MO12212 LB2-01-06194



V1.1